# Autoclavable extension cable

## EN: INSTRUCTIONS FOR USE

BG:	ИНСТРУКЦИИ ЗА УПОТРЕБА	LT:	NAUDOJIMO INSTRUKCIJOS
CS:	NÁVOD K POUŽITÍ	LV:	NORĀDĪJUM PAR LIETOŠANU
DA:	BRUGSANVISNING	MK:	УПАТСТВО ЗА УПОТРЕБА
DE:	GEBRAUCHSANWEISUNG	NL:	GEbruKSAANWIJZING
EL:	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	NO:	BRUKSANVISNING
ES:	INSTRUCCIONES DE USO	PL:	INSTRUKCJA STOSOWANIA
ET:	KASUTUSJUHEND	PT:	INSTRUÇÕES DE UTILIZAÇÃO
FI:	KÄYTTÖOHJEET	RO:	INSTRUCȚIUNI DE UTILIZARE
FR:	CONSIGNES D'UTILISATION	RU:	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
HR:	UPUTE ZA UPORABU	SK:	POKyny NA POUŽITIE
HU:	HASZNÁLATI ÚTMUTATÓ	SL:	NAVODILA ZA UPORABO
ID:	PETUNJUK PENGGUNAAN	SR:	UPUTSTVO ZA UPOTREBU
IS:	NOTKUNARLEIÐBEININGAR	SV:	BRUKSANVISNING
IT:	ISTRUZIONI PER L'USO	TR:	KULLANIM TALIMATLARI
JA:	使用説明書	ZH-S:	使用說明
KO:	사용 설명서	ZH-T:	使用說明書

Technomed Europe  
Amerikalaan 71  
6199 AE Maastricht-Airport  
The Netherlands  
[www.technomed.nl](http://www.technomed.nl)

## EN: Authorized representatives

BG: Упълномощени представители / CS: Autorizovaní zástupci / DA: Autoriserede repræsentanter / DE: Bevollmächtigte Vertreter / EL: Εξουσιοδοτημένοι αντιπρόσωποι / ES: Representantes autorizados / ET: Volitatuud esindajad / FI: Valtuutetut edustajat / FR: Représentants autorisés / HR: Ovlašteni predstavnici / HU: Maghatalmazotti képviselők / ID: Perwakilan resmi / IS: Váðunefndir fulltrúa / IT: Rappresentanti autorizzati / JA: 認定販売代理店 / KO: 공식 대리인 / LT: Įgalioję atstovai / LV: Pilnvarotie pārstāvji / MK: Овластени претставници / NL: Bevoegde vertegenwoordigers / NO: Autoriserte representanter / PL: Autoryzowani przedstawiciele / PT: Representantes autorizados / RO: Reprezentanți autorizați / RU: Уполномоченные представители / SK: Autorizovaní zástupcovia / SL: Pooblašteni zastopniki / SR: Ovlašćeni predstavnici / SV: Auktoriserade representanter / TR: Yetkili temsilciler / ZH-S: 授权代表 / ZH-T: 授權代表

United Kingdom  
The Magstim Company Ltd.  
Spring Gardens  
Whitland SA34 0HR  
United Kingdom

Switzerland  
MedEnvoy Switzerland  
Gothardstrasse 28  
6302 Zug  
Switzerland

Marketed by

**Biosense Webster**

PART OF THE **JOHNSON & JOHNSON** FAMILY OF COMPANIES

# Instructions for use - EN

## AUTOCLAVABLE EXTENSION CABLE

Federal (USA) law restricts this device to sale by or on the order of a physician.

### Indications / Intended use

The cable is designed to connect a Biosense Webster electrophysiology catheter to the appropriate stimulation and / or recording equipment. The cable may also be used as an extension cable to permit catheter connection to equipment that is located out of the immediate area of use. The cable may be reused subject to the reuse instructions provided below.

### Contraindications

There are no known contraindications to this cable.

### Patient target group

These cables are targeted at patients who have been diagnosed with cardiac arrhythmias and are undergoing an electrophysiology procedure.

### Intended users

The intended users of these cables are appropriately trained personnel in a fully equipped electrophysiology laboratory.

### Clinical benefit and product compatibility

These cables provide a means to connect a Biosense Webster electrophysiology catheter to the appropriate equipment. These cables contain 4, 6, 8 or 10 connectors. Use the cable that has the same number of connectors as the diagnostic catheter. For a catheter with 20 connectors, use two (2) cables containing 10 connectors each. The product on its own does not have a clinical benefit and needs to be used in combination with its compatible products to have a clinical benefit.

### Warnings and precautions

- \* Use the cable prior to the "Use By" date on the package label.
- \* The device is delivered sterile using ethylene oxide gas. Some patients may be allergic to residuals of this gas.
- \* In the event the sterile barrier system is opened or damaged prior to use, the sterile state of the product can not be guaranteed. In these circumstances, the device needs to be resterilized prior to use according to instructions.
- \* Patient or operator injury may result from improper handling, improper connection, or improper use of the cable.
- \* If the cable is used in the presence of electrical equipment, noise may be induced into the cable.
- \* If conventional line powered equipment is used near the patient, follow electromechanical device and hospital standards to prevent extraneous electrical currents from reaching the electrophysiological equipment, catheter and heart. Such currents may cause lethal arrhythmias.
- \* Do not use the device in strong electrical fields, as any induced electrical field may influence the stimulation output on the connected equipment (if applicable), may make the readings on the connected equipment unreliable and may result in localized tissue heating.
- \* Do not excessively bend or fold the cable as this may damage the cable. The cable should be visually inspected for any damage (i.e., kinks, nicks, crushed or elongated sections) prior to use or reuse. Do not attempt to repair any damages. In case of doubt, discard the cable and do not reuse.
- \* If a break occurs in the cable wire or if the cable becomes otherwise electrically discontinuous, arcing may occur in the patient-return or active circuit and may burn the patient or create a fire.
- \* Ensure proper cable connection by referring to the connector designations at each end of the cable.
- \* To prevent damage, store the cable in a clean and secure area.

### MRI safety information

Do not use in an MRI scanning environment. The device has not been evaluated for safety in an MR environment. It has not been tested for heating or unwanted movement in an MR environment. Its safety in the MR environment is unknown. Performing an MR exam on a person who has this medical device inserted or positioned on them may result in injury and / or device malfunction.

### Product reuse and sterilization

The cable is packaged sterile from Technomed. The cable may be re-processed and sterilized using steam autoclave or ethylene oxide (EtO) gas. After use, clean the cable thoroughly as described below. The cable can withstand cleaning with soap solutions or alcohol and may be left in contact with cold sterilization solutions for up to two hours.

### CLEANING AND STERILIZATION

The cables are provided sterile. They may be sterilized as described below.

#### Warnings and precautions:

1. Do not immerse the cable connectors in liquid. Do not wipe the cable connectors with liquid. The presence of liquid inside the connectors may result in adverse effects, including improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury.
2. Do not expose the cable to organic solvents because they may damage the cable.
3. Failure to properly clean the cable could lead to inadequate resterilization.
4. When sterilizing the cables, do not exceed 140°C because temperatures higher than 140°C may damage the cable.
5. Automated cleaning of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury.

#### Limitations on cleaning and sterilization:

Clean and sterilize the cable after each use. The maximum number of sterilization cycles is specified in table 1. Exceeding the specified number of sterilization cycles may impair the proper functioning of the cable.

#### Initial treatment at point of use:

There are no requirements for initial treatment at point of use.

#### Containment and transportation:

Place the cable in a clean biohazard bag after use for subsequent cleaning, disinfection and sterilization. Clean, disinfect and sterilize the cable as soon as it is reasonably practical after use.

#### Preparation before cleaning:

Visually inspect the cable for any damage such as crushed or elongated sections, cuts, kinks or nicks. If damage is observed, do not use the cable and dispose of it per the "Disposal" section below.

#### Automated cleaning:

Automated cleaning of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury.

#### Manual cleaning:

1. Wear sterile gloves while handling the cable.
2. Using a wipe\*, use a twisting motion to wipe the entire surface of the cable, but not the cable connectors, to remove soil. Pay particular attention to complex cable features (such as non-smooth areas, joints and crevices).
3. Using another wipe\*\*, use a twisting motion to thoroughly wipe the entire surface of the cable to remove soil.
4. Visually inspect the cable to ensure there is no visible soil. Additional wipes may be used to remove any visible soil.
5. Leave the cable in ambient conditions until it is visibly dry.
6. When the cable is dry, visually inspect the entire surface of the cable in adequate lighting to ensure complete removal of visible soil.

#### \* Recommended wipes:

- A. CaviWipes™ or any other wipes with 17.2% isopropanol and 0.28% disubutylphenoxyethoxyethylmethylbenzyl ammonium chloride may be used.
- B. Clorox Healthcare® Bleach Germicidal Wipes or 70% isopropyl alcohol (IPA) wipes

#### Disinfection acceptable methods:

##### A. CaviWipes™

1. Wear sterile gloves while handling the cable.
2. Ensure manual cleaning (see the "Manual Cleaning" instructions above) is performed before disinfection.
3. Using a wipe\*, use a twisting motion to thoroughly wipe the entire surface of the cable. Pay particular attention to complex cable features (such as non-smooth areas, joints and crevices).
4. Allow the surface of the cable to remain visibly wet for 3 minutes at room temperature (approximately 68°F / 20°C).
5. Additional wipes may be used to ensure that the cable remains visibly wet for 3 minutes.
6. Leave the cable in ambient conditions until it is visibly dry.
7. Visually inspect the entire surface of the cable in adequate lighting to ensure the cable is visibly clean.

\* Recommended wipes: CaviWipes™ or any other wipes with 17.2% isopropanol and 0.28% disubutylphenoxyethoxyethylmethylbenzyl ammonium chloride may be used.

##### B. Clorox Healthcare® Bleach Germicidal Wipes or 70% isopropyl alcohol (IPA) wipes

1. Wear sterile gloves while handling the cable.
2. Ensure manual cleaning (see the "Cleaning: Manual" instructions in this table) is performed before disinfection.
3. Using two (2) wipes\*\*, use a twisting motion to thoroughly wipe the entire surface of the cable. Pay particular attention to complex cable features (such as non smooth areas, joints, and crevices).
4. Using a wipe\*\*, use a twisting motion to wipe the cable and allow the surface to remain visibly wet for the corresponding contact time\*\*\*.
5. Additional wipes may be used to maintain the contact time.
6. Leave the cable in ambient conditions until it is visibly dry.
7. Visually inspect the entire surface of the cable in adequate lighting to ensure the cable is visibly clean.

\*\*Recommended Wipes and

\*\*\*Contact Time: Clorox Healthcare® Bleach Germicidal Wipes: 3 min or 70% isopropyl alcohol (IPA) wipes: 6 min

#### Drying:

Drying is performed as part of the disinfection process (see the "Disinfection" instructions).

#### Maintenance, inspection and testing:

Before packaging, visually inspect the cable. If visible soil is present, repeat the "Manual Cleaning" and "Disinfection" instructions. If there is evidence of damage, do not use the cable and dispose of it per the "Disposal" section below.

#### Packaging:

If the cable is stored after cleaning and disinfection, place it in a breathable pouch that is approved by the FDA and / or by local authorities. Coil the cable into a loop when placing it into the pouch. Do not fold the cable because doing so may break the cable. Use a pouch large enough to prevent stress to the pouch seams and to prevent excessive bending of the cable. A pouch size of 19cm x 33cm (7.5" x 13") or larger is recommended. For storage conditions, refer to the "Storage" section below.

#### Resterilization:

See Table 1 for resterilization methods and parameters. The cable is designed for a maximum of ten sterilization cycles. Resterilize the cable after each use. If the cable is stored after resterilization, the cable should be packaged prior to the sterilization cycle, per "Packaging" section above.\*

The instructions provided above have been validated by the medical device manufacturer as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing as performed (using equipment, materials, and personnel in the processing facility) achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

Table 1 – Re-sterilization parameters for cables

	Exposure Temperature	Exposure Time	Minimum Drying Times	Gas Concentration	Humidity (RH)	Accumulation / Aeration Time at 55°C (131°F)
Gravity Steam Re-sterilization	121°C (250°F)	30 min	15 - 30 min	---	---	---
	132°C (270°F)	15 min	15 - 30 min	---	---	---
	135°C (275°F)	10 min	30 min	---	---	---
Dynamic Air Removal (Pre-Vac) Re-sterilization	134°C (274°F)	3 min	16 min	---	---	---
	135°C (275°F)	3 min	16 min	---	---	---
100% EtO Gas Re-sterilization	55°C (131°F)	1 hr	---	725 mg/L	30% - 80%	12 hr
<p><b>Do not reuse the cable after ten (10) uses.</b> The cable is designed for a maximum of ten sterilization cycles. Exceeding ten uses may impair the proper functioning of the cable. The cable should be discarded after the ten re-sterilization cycles have been reached. Technomed Europe will not be responsible for any direct or consequential damages or expenses that result from reusing the cable beyond the indicated reuse limitations.</p>						

#### STORAGE

Refer to the product label for recommended storage conditions.

#### DISPOSAL

Recycle components or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

#### REPORTING INCIDENTS

Per regulation (EU)2017/745 on medical devices, if this device is used in the European Union (or in a country with an identical regime) and a serious incident occurs during the use of this device or as a result of use of this device, report the incident to the manufacturer Technomed Europe at [quality@technomed.nl](mailto:quality@technomed.nl) and to your national authority.

#### ADVERSE REACTIONS

A number of adverse reactions have been documented for electrophysiology procedures, including heart block, pulmonary vein stenosis, esophageal fistula and / or injury, stroke (cerebrovascular accident), other arrhythmias (outside diagnosis), life-threatening arrhythmias, myocardial infarction, cardiac perforation, pericardial effusion, cardiac tamponade, thrombosis, embolism, pulmonary embolism, air embolism, valvular damage, phrenic nerve injury, vagal nerve injury, pericarditis, coronary artery stenosis, vessel perforation (peripheral and / or central), soft tissue injury, persistent atrial communication, device-related infection, embolization of components, device entrapment, surgical intervention (additional), skin burns, and localized skin reaction.

The following complications associated with cardiac catheterization have also been reported in the literature: major bleed, hematoma, reaction to medications, allergic reaction, vascular access complication, damage to vasculature, implanted device interactions, renal artery stenosis, pneumothorax, ST segment changes, fluid overload, urinary catheter complications, hypotension, sepsis, wound infection, respiratory failure, heart injury, renal injury, heart failure, cardiac arrest and death.

#### ELECTRONIC INSTRUCTIONS FOR USE

This document is available at [www.global-ifu.com/634192](http://www.global-ifu.com/634192)

#### Disclaimer of warranty and limitation of remedy

**There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Technomed Europe products described in this publication. Under no circumstances shall Technomed Europe be liable for any direct, incidental or consequential damages other than as expressly provided by specific law. No person has the authority to bind Technomed Europe to any representation or warranty except as specifically set forth herein.**

Descriptions or specifications in Technomed Europe printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

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