# Probes

**(Disposable stimulation probe (insulated electrode) for neuro monitoring)**

- **Monopolar (Tapered/Flush tip/XL/Central blade nerve)**
- **Bipolar (precision)**
- **Concentric**
- **Hooked/ Curved (double, triple)**
- **Ball tip/Pedicle screw (XL)**

## PROBES

### (Disposable stimulation probe (insulated electrode) for neuro monitoring)

Monopolar (Tapered/Flush tip/XL/Central blade nerve), Bipolar (precision), Concentric, Hooked/ Curved (double, triple), Ball tip/Pedicle screw (XL)

#### Intended use and indications for use

Disposable probes are intended to locate and identify cranial motor nerves during ear-nose-throat (ENT) and intracranial procedures; to locate and identify cranial and peripheral motor nerves during surgery (including spinal nerve roots); and to locate, identify and monitor cranial motor nerves during surgery. The motor nerves are monitored by detecting electromyographical activity in the muscles they innervate.

All of these products are applied during intraoperative neurophysiological monitoring in surgical procedures to locate and identify peripheral, cranial and spinal nerve roots that are at risk of damage. The probe design appropriate for the surgical procedure is determined by the physician.

#### Intended users

The intended users are healthcare professionals specifically trained and certified in electrophysiology techniques.

#### Patient target groups

The patient target groups are any patients identified by these medical experts to benefit from the diagnostic procedures as described above.

#### WARNING

Careless use of any element in electrosurgical systems may cause serious burns. The use of the device with safely functioning and compatible electrosurgical generators and accessories is the user's responsibility. Read and understand all warnings, precautions, cautions and instructions before attempting to use any electrosurgical system.

When stimulating nerves, an appropriate voltage / current should be applied according to the stimulator manufacturer's operating instructions.

#### Cautions

\* Federal law (USA) restricts this device to sale by or on the order of a physician.

\* Do not use the device if the sterile packaging is opened or damaged or if the use-before-date has passed.

\* Inspect the device after opening in a sterile field. Discard if damaged.

\* For single-patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilizing the device can affect its safety, performance and effectiveness and expose patients and users to unnecessary risks such as infection and transmissible disease.

\* Sterility is guaranteed up to expiration date unless package is opened or damaged.

\* The device has been sterilized using ethylene oxide gas. Some patients may be allergic to residuals of this gas.

\* Do not use the device in strong electrical fields, as any induced electrical field may influence the stimulation output on the connected equipment (if applicable), may make the readings on the connected equipment unreliable and may result in localized tissue heating.

\* To avoid damage to the nerve monitor or tissue, the device must be kept separate from active electrosurgical devices and the patient when not in use. A separate holder is recommended.

\* Keep voltage / current settings at lowest possible level to achieve the desired effects.

\* Hold the shaft of the device during bending to avoid damage to the device. If bent more than 30°, the device may break or insulation may be damaged.

\* Do not cut or modify the shaft of the device before or during surgery.

\* Do not stimulate for a prolonged period of time to prevent nerve fatigue.

#### MRI safety information

Do not use in an MRI scanning environment. The device has not been evaluated for safety in an MR environment. It has not been tested for heating or unwanted movement in an MR environment. Its safety in the MR environment is unknown. Performing an MR exam on a person who has this medical device inserted or positioned on them may result in injury and / or device malfunction.

#### Storage, transport and handling

Keep away from extreme temperature, humidity and direct sunlight.

#### Compatibility with accompanying devices

This device is equipped with DIN 42802 connectors and is compatible with any intraoperative neuromonitoring equipment that is equipped with DIN 42802 connections.

#### Disposal

Always discard used devices in a properly marked medical biohazard container.

#### Reporting incidents

All serious incidents associated with the use of this device should be reported to Technomed Europe at [quality@technomed.nl](mailto:quality@technomed.nl) and to a competent authority of the country where the user is established.

#### Electronic instructions for use

Access the most recent electronic instructions for use through the link which is provided on the label of the device.

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Pty Ltd
The Magstim Company Ltd.