# Surgical instruments

## Disposable stimulating:

- RING DISSECTOR
- CURVED NEEDLE
- ELEVATOR
- DISSECTORS

STERILE EO

### EN: INSTRUCTIONS FOR USE

BG:	ИНСТРУКЦИИ ЗА УПОТРЕБА
CS:	NÁVOD K POUŽITÍ
DA:	BRUGSANVISNING
DE:	GEBRAUCHSANWEISUNG
EL:	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES:	INSTRUCCIONES DE USO
ET:	KASUTUSJUHEND
FI:	KÄYTTÖOHJEET
FR:	CONSIGNES D'UTILISATION
HR:	UPUTE ZA UPORABU
HU:	HASZNÁLATI ÚTMUTATÓ
IS:	NOTKUNARLEIÐBEININGAR
IT:	ISTRUZIONI PER L'USO
JA:	使用説明書
KO:	사용 설명서
LT:	NAUDOJIMO INSTRUKCIJOS
LV:	NORĀDĪJUMI PAR LIETOŠANU
MK:	УПАТСТВО ЗА УПОТРЕБА
NL:	GEBRUIKSAANWIJZING
NO:	BRUKSANVISNING
PL:	INSTRUKCJA STOSOWANIA
PT:	INSTRUÇÕES DE UTILIZAÇÃO
RO:	INSTRUCȚIUNI DE UTILIZARE
RU:	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
SK:	POKYNY NA POUŽITIE
SL:	NAVODILA ZA UPORABO
SR:	UPUTSTVO ZA UPOTREBU
SV:	BRUKSANVISNING
TR:	KULLANIM TALİMATLARI
ZH-S:	使用说明

**DISPOSABLE STIMULATING:  
RING DISSECTOR, CURVED NEEDLE, ELEVATOR, DISSECTORS****Intended use and indications for use**

The disposable stimulating ring dissector, stimulating curved needle, stimulating elevator and stimulating dissectors are indicated for tissue dissection and stimulation of cranial and peripheral motor nerves for location and identification purposes during surgery, including spinal nerve roots. They are used in electromyography (EMG).

These devices are applied during intraoperative nerve monitoring in surgical procedures to locate and identify peripheral, cranial and spinal nerve roots that are at risk of damage.

**Intended users**

The intended users are healthcare professionals specifically trained and certified in electrophysiology techniques.

**Patient target groups**

The patient target groups are any patients identified by these medical experts to benefit from the diagnostic procedures as described above.

**Use of product**

\* Remove and discard the protection sheath from the device shaft / tip before use, if applicable.

\* The stimulating ring dissector, stimulating curved needle and stimulating elevator are monopolar probes and require a separate single needle electrode to act as the return. This single needle electrode is supplied in the carton (depending on the variant). The instruments have different tip shapes for performing different dissecting activities.

\* The stimulating dissector is a monopolar probe and requires a separate single needle electrode to act as the return. This single needle electrode is supplied in the carton (depending on the variant). Both tips of the dissector are active. To use the dissector as a probe, close the tips and touch the nerve. If the nerve monitor responds when the tips are apart during dissection, either of the tips may be close to the nerve.

**WARNING**

Careless use of any element in electrosurgical systems may cause serious burns. The use of the device with safely functioning and compatible electrosurgical generators and accessories is the user's responsibility. Read and understand all warnings, precautions, cautions and instructions before attempting to use any electrosurgical system.

When stimulating nerves, an appropriate voltage / current should be applied according to the stimulator manufacturer's operating instructions.

**Cautions**

\* Federal law (USA) restricts this device to sale by or on the order of a physician.

\* Do not use the device if the sterile packaging is opened or damaged or if the use-before-date has passed.

\* Inspect the device after opening in a sterile field. Discard if damaged.

\* For single-patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilizing the device can affect its safety, performance and effectiveness and expose patients and users to unnecessary risks such as infection and transmissible disease.

\* Sterility is guaranteed up to expiration date unless package is opened or damaged.

\* The device has been sterilized using ethylene oxide gas. Some patients may be allergic to residuals of this gas.

\* Do not use the device in strong electrical fields, as any induced electrical field may influence the stimulation output on the connected equipment (if applicable), may make the readings on the connected equipment unreliable and may result in localized tissue heating.

\* To avoid damage to the nerve monitor or tissue, the device must be kept separate from active electrosurgical devices and the patient when not in use. A separate holder is recommended.

\* Keep voltage / current settings at the lowest possible level to achieve the desired effects.

\* Do not stimulate for a prolonged period of time to prevent nerve fatigue.

**MRI safety information**

Do not use in an MRI scanning environment. The device has not been evaluated for safety in an MR environment. It has not been tested for heating or unwanted movement in an MR environment. Its safety in the MR environment is unknown. Performing an MR exam on a person who has this medical device inserted or positioned on them may result in injury and / or device malfunction.

**Storage, transport and handling**

Keep away from extreme temperatures, humidity and direct sunlight.

**Compatibility with accompanying devices**

This device is equipped with DIN 42802 connectors and is compatible with any intraoperative neuromonitoring equipment that is equipped with DIN 42802 connections.

**Disposal**

Always discard used devices in a properly marked medical biohazard container.

**Reporting incidents**

All serious incidents associated with the use of this device should be reported to Technomed Europe at [qualitytechnomed@welcony.com](mailto:qualitytechnomed@welcony.com) and to a competent authority of the country where the user is established.

**Electronic instructions for use**

Access the most recent electronic instructions for use through the link which is provided on the label of the device.

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Pty Ltd
The Magstim Company Ltd.